



ELECTRIC HI-LO TREATMENT TABLE

MODEL 4707

Shown in left hand view with options:

06 - Backrest

27 - Paper Cutter & Dispenser



Model 4707

Electric Hi-Lo Treatment Table



- Adjustable height motorized hydraulic cylinders are built into each (4) corner steel leg.
- "Hands Free" foot control activation.
- Apron and shelf are almond laminate.
- Black powder coated legs with levelers.
- 110V, 60Hz, 2.4 Amps.
- Choice of (18) standard vinyl colors.

OPTIONS

- 06* - Backrest, adjustable
- 07* - Nose Cutout
- 27 - Paper Dispenser & Cutter Combo
- 99* - Left Hand Table Model



*Factory Installed

MODEL	LENGTH	WIDTH	HEIGHT
4707	72"	30"	27" - 37"