



PRODUCTS BY SPECIALTY: CARDIOLOGY

05 – Headrest is now standard



Model 4790 with Drop-Away "Closed"
Shown with options 074 – 4" Locking Casters (not intended for patient transport) and 27 – Paper Dispenser & Cutter



A. Model 4790 with Drop Away "Open"

COMMON FEATURES

ECHO-SCAN™ ECHOCARDIOGRAPH TABLES MODELS 4790, 4892



- 3" High-Density urethane foam top.
- 10" wide x 12" deep drop-away section on durable piano hinge.
- Sliding latch activates drop-away section.
- Choice of (18) standard vinyl colors

A. Model 4790 Powermatic® Echo-Scan™ Table

- Power height adjustment from 20" wheelchair height to 29" high.
- "Hands-free" foot control.
- Hospital grade plug.
- 110V, 60Hz, 2.8 Amps. U.L. Listed.
- Includes manually adjustable headrest.

MODEL	LENGTH	WIDTH	HEIGHT
4790	76"	27"	20"- 29"

B. Model 4892 Echo-Scan™ Table

- Two laminated side panel legs.
- Convenient laminate storage shelf.
- Folkstone Gray laminate finish.

MODEL	LENGTH	WIDTH	HEIGHT
4892	76"	27"	31"

OPTIONS

FOR MODELS 4790, 4892

- *04 – Table Safety Strap
- *05 – Headrest, adjustable (4892 only)
- 27 – Paper Dispenser & Cutter Combo
- 36 – Wedge Pillow
- *074 – (4) 4" Locking Casters (4790 only) (adds 4" to unit height)

*FACTORY INSTALLED



B. Model 4892 with Drop-Away "Open"

CONTACT: →

EXAM TABLES

4423

PROFESSIONAL STOOLS

2153

CHROMEWARE

2182

2188

BLOOD CHAIRS

4423

PRO-LINET™ CABINETRY

B-401

TREATMENT TABLES

4024

FOOT STOOLS

2038