



DOC-IN-A-BOX



D. Model 2182



A. Model 4002-027
Shown with Option 27 Paper Dispenser & Cutter Combo



B. Model 9038-20-346



C. Model 2030



A. 4002-027 Treatment Table H-Brace support is recessed into table leg with unique 4-sided "lock-tite" joint. High pressure oak laminate legs, apron, and stretchers. 2" High-density foam top. Resist-All upholstery finish. Latex free. Paper dispenser and cutter are optional. 400 lbs. capacity.

B. 9038-20-346 Mobile Cart Maintenance-free laminated surfaces on all surfaces. High impact polystyrene drawer interiors with rounded corners. Oak laminate cabinets with satin aluminum wire drawer pulls. 2 1/2" twin wheel hooded swivel nylon casters. 100 lbs. weight capacity. 18 1/2"W x 18 1/2"D x 30"H.

C. 2030 Footstool with Safety Handrail 34"H Chrome handrail. Non-slip treaded step 14" x 11" x 9"H. Chrome steel tubular lags with safety tips.

D. 2182 Gooseneck Lamp 12"L flexible gooseneck arm on 12"L adjustable telescoping tube. Inner and outer shade dissipates heat. Heavy plated steel base.

More Quality Products from Hausmann: Full product details can be found in the Hausmann catalog or online at www.hausmann.com



4416 600 LBS



4423
Available options include: Safety Rail, Pelvic Tilt, and Irrigation Pan



8728-346



9018-10-346



2000



7031
Shown with optional 05 Headrest and 27 Paper Dispenser & Cutter

CONTACT: →