



HAUSMANN MOST POPULAR

TREATMENT TABLES



A. Model 4002



A. Series 4002

- H-brace is recessed into table leg with unique 4-sided "Lock-Tite" joint.
- High pressure laminate legs, apron and stretchers.



B. Model 4024



B. Series 4024 w/Shelf

- H-brace with laminate storage shelf.
- High pressure laminate legs, apron and stretchers.

 = Ships out 3 days from order.*



C. Model 4041



C. Series 4041 w/Storage

- H-brace with laminate shelf and open-shelf storage cabinet.
- High pressure laminate legs, apron and stretchers.



D. Model 4043



D. Series 4043

- Fully-enclosed, all-laminate storage cabinet with sliding laminate doors.
- Interior has divider and one adjustable shelf.
- Black laminated plywood base.

Common Features

- Comfortable 2" thick high-density foam upholstered tops.
- Laminate: Natural Oak, Wild Cherry or Folkstone Gray
- Upholstery: 18 standard vinyl colors
- Sizes: 72"x24", 72"x27", 72"x30", 78"x30"

E. Model 4075



Wall Folding Treatment/Changing Table

- Sturdy, floor-supported wall mounting plate.
- Two folding high pressure laminate legs.
- Projects only 8" when folded.
- Natural Oak laminate only.



E. Model 4075

Model 4075 folded

*See website for more information.

TREATMENT/EXAM TABLES



A. Model 4795

A. Model 4795 Bariatric Hi-Lo Treatment Table with Power Backrest

- Extra-wide 32" top.
- Power height elevation from 21" wheelchair height to 30" high.
- Power backrest elevates from 0° to 75°.
- Convenient hand control operation.
- Black laminate frame and base plate with aluminum trim and levelers.
- Upholstery: 18 standard vinyl colors.



C. Model 4416

C. Model 4416 X-L Power-Back Bariatric Exam Table

- Extra-wide 30" top.
- Backrest is motorized 0° to 75°, with hand control.
- Spacious footstep with anti-slip treaded surface and front supports.
- Concealed stirrups.
- Slide-out 20"W legrest with pad extends 18".
- Pelvic tilt with dual side controls.
- Paper dispenser and cutter.
- Upholstery: 18 standard vinyl colors.



B. Model 4750

B. Model 4750, 4756, 4757 Powermatic® Table

- The Powermatic® direct drive mechanism provides an infinite number of height positions.
- Patients can mount or dismount without a footstep.
- Power backrest elevates from 0° to 75°.
- 2" high-density urethane foam vinyl top.
- "Hands free" foot control operation.
- Includes paper dispenser and cutter.
- Upholstery: 18 standard vinyl colors.
- Model 4750 with Flat Top
- Model 4756 with Gas-Spring Backrest
- Model 4757 with Power Backrest and foot control.



D. Model 4746

D. Model 4746, 4747 Powermatic® Treatment Table

- Power height adjusts from 26" to 38".
- Gas-spring pneumatic 30" backrest elevates to 75° with lever control.
- Hands-free foot control operation.
- Laminate: Folkstone Gray or choice of any high pressure laminate solid color or woodgrain.
- Upholstery: 18 standard vinyl colors.
- Model 4746 with Gas-Spring Backrest.
- Model 4747 with Power Backrest and foot control.





A. Model 4723

A. Model 4723 ADA Hi-Lo Power Plinth Table

- 72"L x 30"W
- Motorized table lowers to 20" for wheelchair transfer. Elevates up to 30" high.
- Electrically activated hydraulic cylinders are built into corner telescoping legs.
- Hands-free foot control operation.
- Black laminate storage shelf.
- Upholstery: 10 Pro-Form™ vinyl colors.



B. Model 7604-752

B. Model 7604 Portable Massage Treatment Table

- 73"L x 28"W
- 2½" thick multi-layer semi-firm foam top.
- Push button height adjustment from 24"-32".
- Precision welded, light weight aluminum frame weighs only 33 lb.
- Includes adjustable face cradle with ultra soft face cushion, and carry bag.
- Comfortable PVC-free heavy-weight urethane upholstery: 752 Blue.



C. Model 4945



C. Model 4945 Pediatric Table

- Wood has no added urea-formaldehyde for safer indoor quality. PVC-free construction.
- Digital scale with lb-oz or kg readout. Reads in 1 oz increments.
- 39"L infantometer with 10"L sliding fold-away flap.
- Cabinet has two drawers and two doors.
- Laminate: Folkstone Gray or choice of any high pressure laminate solid color or woodgrain.
- Upholstery: 5 PVC-free urethane colors.



D. Model 2151

400
LBS

E. Model 2153

250
LBS



F. Model 2157

250
LBS



G. Model 3415

65
LBS

Common Features

MODELS 2153, 2157

3

- 16" diameter seat, 3½" thick foam.
- Twin wheel nylon casters.
- Upholstery: 18 standard vinyl colors.

D. Model 2151 Heavy Duty Air-Lift Stool

- Heavy-duty pneumatic adjustment with D-handle.

E. Model 2153 Economy Air-Lift Stool

- Single lever height control.

F. Model 2157 Economy Air-Lift Stool

- 360° circular ring control handle.

G. Model 3415 Easy-Lift Over Bed Table

- Ultra-smooth gas lift adjustment, locks table at desired height.

TILT TABLES/MAT PLATFORMS

A. Model 6045



A. Model 6045 Bariatric Electric Tilt Table



- Extra-wide 34" x 77" long top upholstered in 709 Gray vinyl with 2" thick high-density urethane foam.
- 23" extra-wide x 18" deep fixed foot board.
- Three extra-long patient safety straps.
- Tilt indicator displays angle of tilt 0° to 90°.
- Fixed table height of 33".
- Hand control operable from either side of the table.
- Not designed for cardiac testing use.


B. Model 6058



B. Model 6058 Wheelchair Accessible Hi-Lo Tilt Table



- Electric Hi-Lo lift mechanism adjusts from low 19" wheelchair height to 33" high.
- Tilts from horizontal to 90°.
- Hands-free foot control operation.
- Includes emergency tilt release handle.
- Four 3" diameter swiveling locking casters.
- Heavy-duty steel frame is a cream powder coated epoxy finish.
- Not designed for cardiac testing use.



C. Model 1440-57-113

C. Series 1440 Dual Lift Powermatic® Mat Platforms



- Power height adjustment from 19" wheelchair height to 28" high.
- Two direct drive, synchronized lift mechanisms for unsurpassed power.
- Tear resistant, seamless nylon-reinforced vinyl upholstered top with 2" thick high-density urethane foam padding for comfort and rounded corners for safety.
- Hands free foot control.
- Available in four sizes and four vinyl colors.



D. Model 1427-57-121

D. Series 1427 Bariatric Mat Platforms



- Motorized lift system can securely lift 750 lb from 20" wheelchair height to 30" high.
- Steel frame and 4 corner legs have durable black powder coat finish.
- 4-leg "open" design is patient sling-lift accessible on all four sides.
- Tear resistant, seamless nylon-reinforced vinyl upholstered top with 2" thick high-density urethane foam top.
- Hands-free foot control.
- Available in four sizes and four vinyl colors.

A. Model 6175



A. Model 6175

Hi-Lo Econo-Line Stand-In Table with Electric Patient Lift

- Motorized patient lift from sitting to standing position.
- Motorized height adjustment of durable Sail White laminate top with cutout.
- Includes lift support harness.
- Hand control operation.
- Padded knee, chest and back supports.
- Cream powder coated heavy gauge steel frame.
- Locking 3" diameter heavy-duty casters.

B. Model 6185



B. Model 6185

Bariatric 500 lb Electric Stand-In Table

- Unit power lifts 500 lb patient from sitting to standing position.
- Motorized height adjustment of extra-wide raised rim laminate top with cutout.
- Hand control operation.
- Padded knee, chest and back supports.
- Includes bariatric and standard size patient lift harnesses.



D. Model 4375

D. Model 4375

Electric Work Table



- Electric lift mechanism elevates top.
- Hand control mounts under work surface.
- Ideal height range for sitting, standing or wheelchair use.
- Easily accommodates four ambulatory patients or two wheelchairs and therapist(s).
- Natural Oak laminate top with black T-molding.
- Black laminate base plate with leveling legs.



E. Model 6674

E. Model 6674

Horseshoe Therapy Table



- U-shape top accommodates up to four wheelchair patients and therapist.
- Top is only 18" deep for easy access.
- Durable, Sail White high pressure laminate top.
- Five adjustable height metal legs.
- Hand control knob on each leg allows 1" adjustments from 26" to 34" high.

WEIGHT RACKS



A. Model 5560-100



B. Model 5562-100

A. Model 5560-100 Multi-Purpose Combination Rack with Accessories



- Distortion-free glass mirror with ANSI “shatter-stop” safety backing.
- 3” diameter swiveling casters.
- Includes:
 - (1) Set of (16) Hugger® Weights, (1) Set of (5) REP™ Bands, (2) Sets of (10) dumbbells.

Common Features

- Built-in dispenser for REP™ Bands.
- Pegboard with eleven safety hooks.
- Hooks for dumbbells(20).
- Base plate with bumper T-molding.
- Black or Natural Oak laminate finish.

C. Model 5556-100 Econo Weight Rack with Accessories



- Built-in compartment for weight bars.
- 3” diameter swiveling casters.
- Includes:
 - (1) Set of (16) Hugger® Weights, (2) Sets of (10) dumbbells, and (1) Set of (3) Weight Bars.



C. Model 5556-100

B. Model 5562-100 Mega Rack with Accessories



- Convenient drop lid work counter with lock.
- Three 8” high laminate drawers.
- Built-in dispensers for VELCRO® rolls.
- 3 sets of hooks for Weight Bars.
- Distortion-free glass mirror with ANSI “shatter-stop” safety backing.
- 4” diameter twin wheel, locking casters.
- Includes:
 - (1) Set of (16) Hugger® Weights,
 - (1) Set of (5) REP™ Bands, (2) Sets of (10) dumbbells, and (1) Set of (3) Weight Bars.

D. Model 5592-100 Compact Weight Rack with Accessories



- Storage compartment for weight bars.
- Laminate shelves with 3 plastic storage tubs.
- 3” diameter swiveling casters.
- Includes:
 - (1) Set of (16) Hugger® Weights, (1) Set of (5) REP™ Bands, (2) Sets of (10) dumbbells, (1) Set of (5) Weight Bars, and (1) Set of (6) Soft Ball hand weights.



D. Model 5592-100



A. Model 1580



A. Model 1580 ADA Convertible “Space-Saver” Staircase

- Handrails and uprights are satin finished premium hardwood and shipped ready to assemble.
- Ergonomically designed handrails have 12” A.D.A. overhang.
- Compact footprint takes up less space and fits in conventional 8’ high ceiling.
- BriteStep™ has yellow hi-visibility, hi-impact 3mm edging for enhanced safety & durability.
- Can be installed in straight line or 90° corner arrangement.
- A.D.A. 11” deep steps x 36” wide have “anti-slip” safety treads.
- (3) 6” high steps on one end and (4) 4½” high steps on the other end.

B. Model 1398 Patented Mobility Platform with Electric Height Bars



- Motorized height adjustment, 29” to 44” high.
- (4) Direct-drive synchronized lifting columns.
- Push button controls and digital bar height readout.
- Width adjusts manually from 15” to 28”W.
- Grommeted holes on platform for placement ladder and balance beam.
- Placement ladder consists of (11) variably placed ladder steps.
- Balance beam is equipped with “anti-slip” surface and dowel pins: 6’5”L x 6”W x 3”H.



B. Model 1398

Common Features

MODELS 1390, 1396, 1398

- Satin-finish hardwood platform with tapered hardwood ends for easy wheelchair access.
- Two “anti-slip” treads on each end.
- 1½” diameter one-piece stainless steel handrails.

Model 1396 Power Height Parallel Bars (not shown)



- Same features as Model 1398 except platform has no dowel holes and does not include ladder steps or balance beam.

C. Model 1390, 1391, 1392 Height & Width Adjustable Parallel Bars



- Width adjusts from 15” to 28”.
- Heavy gauge black powder coated steel uprights and fittings.
- Each upright telescopes up in 1½” increments and locks into ten height positions.
- 3mm platform edge for greater durability.
- 1390 - 10 feet
- 1391 - 7 feet
- 1392 - 12 feet



C. Model 1390

PARALLEL BARS/MULTI-PURPOSE TABLES/ STOOL/MIRRORS



A. Model 1357



B. Model 6293
w/ DRAWER & CASTERS

A. Model 1357 Electric Height & Width Parallel Bars




- Individual digital controls for both height and width adjustment.
- Battery powered.
- Bars overhang platform by 12" for wheelchair assistance.
- 600 lb weight capacity, 300 lb each side.
- Easy grip 2" bars.
- Width: 19"- 33" • Length: 10'
- Height: 24"- 40"



B. Model 6292, 6293, 6294 Multi-Purpose Tables



- Dual gas spring columns, combined with a balanced table top allows the top to be effortlessly lifted to required height.
- Dual hand control levers on each side releases and lowers the height to desired level.
- Brushed aluminum dual telescoping legs with black powder coated U-shaped base.
- 1" thick, high pressure laminate tops, finished in Folkstone Gray with black T-molding edge banding and rounded corners.
- Model 6293 with Drawer
- Model 6294 with Cutout



C. Model 2418

C. Model 2418 Heavy Duty Color-Coded Nested Footstools



- 3¾", 6¼", and 8¾" color coded heights.
- 19" x 12" platforms with "anti-slip" tread.
- Light weight construction – ONLY 5 lb each.
- All aluminum welded construction.
- Wide "anti-tip" bases with rubber feet.
- "Outrigger" legs are very stable – Great for geriatric & pediatric use!



D. Model 1671 Portable Single Mirror



- Distortion-free glass mirror with ANSI "shatter-stop" safety backing.
- Oak laminate frame and base.
- Finished oak laminate back.
- Mounted on easy rolling 2" diameter swiveling casters.
- Mirror height provides full body image.



WWW.HAUSMANN.COM

CAT 1018 © 2019 HAUSMANN INDUSTRIES

The product dimensions and specifications in this catalog apply to the specific products and models described at the date of issue. Under our policy of continuous product improvement, all dimensions and specifications are subject to change without notice. Visit www.hausmann.com for the most current information and products available.